

Easy-Loc Mounting System for **ZOLL**, M,E,R, X Series & Propaq Devices





Product Catalogue 2020



Surface Mounting Bracket With Rotate Part Number SES/ELOO1

"The extremely robust Easy-loc bracket "The extremely popular Easy-loc bracket has been designed to offer a quick locking facility when storing the defibrillator along with an easy release mechanism giving the user a fast release to suit any emergency. This durable bracket has been designed and tested to suit the high demands of the EMS environment. The Easy-loc carriage box is designed to mount the Zoll E, M, R or X Series Defibrillators. This bracket allows the user to rotate the device to have the perfect view of the screen anywhere on the vehicle. This surface mounting bracket allows the turntable to lock into a stow position to prevent the bracket moving when the Defibrillator is not in use. Handle colour or shape is optional to suit each service requirements. This comes complete as a kit with a Device Mounting Plate. Style of baseplate to suit dedicated device to be stipulated at the time of order.

Overall Dimensions = 295mm X 295mm X 73mm Max Height (Including carriage box).











Fixed Surface Mounting Bracket

Part Number SES/ELOO2

"The extremely robust Easy-loc bracket "The extremely popular Easy-loc bracket has been designed to offer a quick locking facility when storing the defibrillator along with an easy release mechanism giving the user a fast release to suit any emergency. This durable bracket has been designed and tested to suit the high demands of the EMS environment. The Easy-loc carriage box is designed to mount the Zoll E, M or X Series Defibrillators. The surface mounting bracket allows the Defibrillator to store in a fixed position where rotation of the defibrillator screen is not needed. These are ideal for RRV's or ambulances with a small shelf along the trauma wall. Handle colour or shape is optional to suit each service requirements. This comes complete as a kit with a Device Mounting Plate. Style of baseplate to suit dedicated device to be stipulated at the time of order.

Overall Dimensions = 290mm X 200mm X 55mm Max Height (Including carriage box).









X Series Fixed Wall Mounting Bracket

Part Number SES/XELOO3

"The extremely robust Easy-loc bracket "The extremely popular Easy-loc bracket has been designed to offer a quick locking facility when storing the defibrillator along with an easy release mechanism giving the user a fast release to suit any emergency. This durable bracket has been designed and tested to suit the high demands of the EMS environment. The Easy-loc carriage box is designed to mount the Zoll X Series Defibrillator. The wall mounting bracket allows the Defibrillator to store in a fixed position where rotation of the defibrillator screen is not needed. These are ideal when wall mounting of equipment to utilise the vehicle space. This bracket is firmly fixed to the trauma wall. The fixing holes can be moved on our shelf to accommodate any existing mounting holes where original defibrillator or monitors have been mounted. This will ensure that your vehicle is still CEN compliant. This bracket also houses the power adaptor cable to charge the X Series. Handle colour or shape and cushion colour is optional to suit each service requirements. This comes complete as a kit with a Device Mounting Plate. Style of baseplate to suit dedicated device to be stipulated at the time of order.

Overall Dimensions = 250mm X 280mm X 52mm Max Height (Including carriage box).

Backing plate overall height 300 mm.







Sliding Tray Bracket with Rotate

Part Number SES/ELSTOO1 Tray Version Only Part Number SES/ELSTOO2 Wall Mounted Version

"The extremely robust Easy-loc bracket "The extremely popular Easy-loc bracket has been designed to offer a quick locking facility when storing the defibrillator along with an easy release mechanism giving the user a fast release to suit any emergency. This durable bracket has been designed and tested to suit the high demands of the EMS environment. The Easy-loc carriage box is designed to mount the Zoll E, M and X Series Defibrillators. This product is very versatile and can be mounted virtually anywhere. Recommended mounting on a small shelf area at the side of the vehicle (Tray Version) or if the device is wall mounted at the foot end of the stretcher. This allows the user to extend the bracket 180mm forward over the stretcher and back to a locked stow position, saving the user having to stretch over and mount the defibrillator onto the bracket. This bracket allows the user to rotate the device to have the perfect view of the screen anywhere on the vehicle. This surface mounting bracket allows the turntable to lock into a stow position to prevent the bracket moving when the Defibrillator is not in use. Handle colour or shape is optional to suit each service requirements. Style of baseplate to suit dedicated device to be stipulated at the time of order.

Overall Dimensions = 360mm X 345mm X 90mm Max Height (Including carriage box Tray height).



Wall Mounted version shown with Zoll X Series

Tray Mounted version



E,M & X Series Wall Mounting Bracket

Part Number SES/ELOO3 Fixed

Part Number SES/ELOO4 Rotate

"The extremely robust Easy-loc bracket "The extremely popular Easy-loc bracket has been designed to offer a quick locking facility when storing the defibrillator along with an easy release mechanism giving the user a fast release to suit any emergency. This durable bracket has been designed and tested to suit the high demands of the EMS environment. The Easy-loc carriage box is designed to mount the Zoll E,M & X Series Defibrillator. The wall mounting bracket allows the Defibrillator to store in a fixed position where rotation of the defibrillator screen is not needed. These are ideal when wall mounting of equipment to utilise the vehicle space. This bracket is firmly fixed to the trauma wall. The fixing holes can be moved on our shelf to accommodate any existing mounting holes where original defibrillator or monitors have been mounted. This will ensure that your vehicle is still CEN compliant. This bracket also houses the power adaptor cable to charge the Zoll device. Handle colour or shape and cushion colour is optional to suit each service requirements. This comes complete as a kit with Device Mounting Plate. Style of baseplate to suit dedicated device to be stipulated as the time of order.



Easy-Loc Wall Bracket with Rotate







EASYLOC CARRIAGE BOX ASS'Y ONLY

PART NO SES/ELCBA001

"The extremely robust Easy-loc bracket "The extremely popular Easy-loc bracket has been designed to offer a quick locking facility when storing the defibrillator along with an easy release mechanism giving the user a fast release to suit any emergency. This durable bracket has been designed and tested to suit the high demands of the EMS environment. The Easy-loc carriage box is designed to mount the Zoll E, M or X Series Defibrillators.

Overall Box Dimensions = 238mm X 195mm X 45mm






Baseplates for Devices

M-Series

To complement the robust Easy-loc bracket, we manufacture a lightweight aluminium mounting plate to accommodate the M,E,R and X-Series defibrillator. All come as standard with heavy duty Stainless Steel feet for ease of loading the Zoll device into the designated bracket. The Plate colour is optional. Please ask for details.



PART NO SES/ELMP001

EASYLOC M SERIES PLATE WITHOUT NIBP TO FIT XTREME PACK



PART NO SES/ELMP002

EASYLOC M SERIES PLATE WITH NIBP TO FIT XTREME PACK



E-Series



PART NO SES/ELEP001

E SERIES EASYLOC PLATE TO SUIT SOFT BAG INC NYLON RUNNERS

X-Series & Propag





PART NO SES/ELXP001

EASYLOC X SERIES STANDARD PLATE



X-Series & Propaq Continued





PART NO SES/ELXP002 EASYLOC X SERIES PLATE WITH EXTENDED CENTRES

This allows for an extra bag to be fitted to the back of the X-Series when wall mounted.





PART NO SES/ELXP003

EASYLOC X SERIES PLATE WITH BAG SUPPORT LUGS (FOR RESCUE AND MEDICAL BAGS) & STRENGTHENER ADAPTOR PLATE



X-Series & Propag Continued

Adaptor / Strengthener Plate only



PART NO SES/XSAP001

X SERIES HEAVY DUTY ADAPTOR PLATE




PART NO SES/XSAP002

X SERIES ADAPTOR PLATE EXTENDED CENTRES



Handle Options



Handle style options available for the Easy-Doc range of brackets as pictured above. These can be supplied in a colour to suit your exact requirements. Please ask for further details. Handle style preference to be discussed at the time of purchase.











We can manufacture a bespoke bracket solution tailored to suit the exact needs of each ambulance service from fixed, to wall mount with an optional rotate system. Please enquire with your exact needs and let us find the perfect solution. We have vast amounts of experience working alongside many ambulance services and ambulance builders worldwide. We have been proudly supporting Zoll with our mounting solutions since 2003.

For further information and pricing on our mounting solutions please contact us.

Shaun Watkin (Director)

Stafford
Engineering
IN PRECISION, PRODUCTION AND PROTOTYPE PROJECTS
STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS STATEMENT OF THE PROJECTS S

Gainsborough Works St Patrick's Place Stafford ST16 2PN

TEL: +44 (0) 1785 220900 FAX: +44 (0) 1785 244897 MOB: +44 (0) 7812056718

EMAIL: sales@staffordengineering.co.uk or shaun@staffordengineering.co.uk

www.staffordengineering.co.uk





